



Rüsch®  
**TracheoFix™**  
 Curved Tracheostomy Tube Sets

**Supports Effective Longer-Term Tracheostomy Care**

The Rüsch® TracheoFix™ Tracheostomy Tube portfolio from Teleflex is ideally suited for ‘second tube in’ situations, or for mid- to longer-term acute care, rehabilitation, or home care, and is designed to meet the changing needs of a tracheostomy patient. The portfolio includes single-use or single-patient-use kit configurations with cuffed and uncuffed tracheostomy tubes available with or without multiple fenestrations. Moreover, a large range of accessories is available both with sets and separately. The Rüsch® TracheoFix™ Tracheostomy Tube sets are packaged in a handy pouch for bedside use to ensure they are easily accessible in a ward or home care environment, and to protect accessories and the replacement tube from airborne contaminants.

**Designed with Patient Safety and Carers in Mind**

The Rüsch® TracheoFix™ Tracheostomy Tubes have size-specific colour-coded inner cannulas (available with or without a 15 mm connector) that correspond to the matching coloured information on the flange. Moreover, the inner cannula can be removed and cleaned or replaced, helping to prevent the accumulation of biofilm and secretions, which could occlude the outer cannula and restrict ventilation. The flexibility of the flange enables it to conform to the anatomy and provides ease of access to the stoma site for cleaning and dressing changes, while the curved tube is designed to be easy to insert, supporting carers in the community setting.

**Benefits**



**The Clinician**

Broad product portfolio ideally suited for mid- to longer-term acute care, rehabilitation, or home care



**Your Institution**

Single-patient-use devices provide a cost-efficient\* option for tracheostomy



**The Patient**

Design features support the changing needs of patients transitioning from hospital to home

\*Cost efficiencies may vary based on a number of factors, including specific markets, services, and rates.

## Rüsch® TracheoFix™ Curved Tracheostomy Tubes Fast Facts

Accessories and replacement tubes stored and protected in a handy pouch.

Thin-walled tube supports reduced work of breathing and a minimum OD.

Multifunctional accessories and sets with rotary connection locking mechanism – kits dependent on tube selection.

Compatible with the Rüsch® PercuTwist™ Percutaneous Tracheostomy Dilator.

### Rüsch® TracheoFix™ Tracheostomy Tube Sets

	TFCP SET	TFP SET	TFC SET	TF SET
<b>ITEM NUMBER</b>				
Acute care	121410	121420	121400	121415
Home care	859001	859003	859000	859002
<b>FEATURES</b>				
Quantity	5	5	5	5
Order size in mm	7.0 – 11.0	7.0 – 11.0	7.0 – 11.0	7.0 – 11.0
Cuffed tracheostomy tube	✓	–	✓	–
Un-cuffed tracheostomy tube	–	✓	–	✓
Fenestrations	✓	✓	–	–
Inner cannula, 15 mm connector, pre-mounted fixation ring (long)	✓	✓	set of two	✓
Inner cannula, low-profile connector (short)	✓	✓	–	✓
Inner cannula, 15 mm connector, phonation, pre-mounted fixation ring, (long)	✓	✓	–	–
Inner cannula, low-profile connector, speaking valve, (short)	✓	✓	–	–
Fixation ring	✓	✓	–	✓
Cough cap	✓	✓	–	✓
Speaking valve with O <sub>2</sub> adaptor	✓	✓	–	–
Comfort neck band	✓	✓	✓	✓
Insertion aid (obturator)	✓	✓	✓	✓

TFPC, Rüsch® TracheoFix™ Cuffed Phonations (Fenestrations); TFP, Rüsch® TracheoFix™, Uncuffed Phonations (Fenestrations); TFC, Rüsch® TracheoFix™ Cuffed; TF, Rüsch® TracheoFix™, Uncuffed

### Rüsch® TracheoFix™ Tracheostomy Tube Dimensions

ORDER SIZE	ID (MM)	OD (MM)	LENGTH (MM)	CURVE APPROX.	CUFF Ø (MM)
7	7.2	9.0	70.0	78°	27.0
8.5	8.7	10.3	78.0	86°	27.0
10.0	10.2	12.3	86.0	90°	32.0
11.0	11.2	13.3	86.0	88°	35.0

### Jackson Size Conversion Chart

JACKSON			RÜSCH® TRACHEOFIX™		
SIZE	ID (MM)	OD (MM)	SIZE	ID (MM)	OD (MM)
4.0	5.0	9.4	7.0	5.5	9.0
6.0	6.4	10.8	8.5	7.0	10.3
8.0	7.6	12.2	10.0	8.5	12.3
10.0	8.9	13.8	11.0	9.5	13.3

Rüsch® TracheoFix™ ID is shown with inner cannula in situ.

### Rüsch® TracheoFix™ Inner Cannula

ITEM NUMBER	DESCRIPTION	ORDER SIZE	ID (MM)	OD (MM)	COLOUR CODE	CASE QTY
121425	15 mm connector, premounted fixation ring	7.0	5.5	6.7	Blue	10
		8.5	7.0	8.2	Green	
		10.0	8.5	9.7	Red	
		11.0	9.5	10.7	Yellow	
121430	Low profile connector, fixation ring	7.0	5.5	6.7	Blue	10
		8.5	7.0	8.2	Green	
		10.0	8.5	9.7	Red	
		11.0	9.5	10.7	Yellow	
121435	15 mm connector, premounted fixation ring, phonation	7.0	5.5	6.7	Blue	10
		8.5	7.0	8.2	Green	
		10.0	8.5	9.7	Red	
		11.0	9.5	10.7	Yellow	
121440	Low profile connector, fixation ring, phonation	7.0	5.5	6.7	Blue	10
		8.5	7.0	8.2	Green	
		10.0	8.5	9.7	Red	
		11.0	9.5	10.7	Yellow	

### Rüsch® TracheoFix™ Accessories

ITEM NUMBER	DESCRIPTION	ORDER SIZE	CASE QTY
121475	Insertion aid	7.0	1
		8.5	1
		10.0	1
		11.0	1
121445	Sealing cap for weaning	–	4
121450	Phonation valve for 15 mm connector with SoftMembrane technology and with O <sub>2</sub> adaptor	–	2
121455	Phonation valve with SoftMembrane technology and with O <sub>2</sub> adaptor	–	2
121460	Cough cap	–	4
121465	Fixation ring to secure the Rüsch® TracheoFix™ Inner Cannula to low-profile connector	–	4



Contents of unopened, undamaged package are clean, ready to use and for single-patient use. Refer to package insert for current warnings, indications, contraindications, precautions, and Instructions for Use.

Teleflex, the Teleflex logo, Rüsch, PercuTwist, and TracheoFix are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current at time of printing (06/2018). Subject to technical changes without further notice. © 2018 Teleflex Incorporated. All rights reserved. MC-XXXXXX Rev.1