

# ACT2™ Container System

Sterile integrity, combined with handling convenience and long-term economy

**Pilling®**

The ACT2™ by Pilling is the total solution sterilization container system combining the finest materials, design and craftsmanship.

## **LIGHTWEIGHT, STRONG**

Crafted from fade- and scratch-resistant, clear hard-coated aluminum, ACT2™ Sterilization Container Systems are lighter, more space efficient and stronger than other systems.

## **STERILE INTEGRITY**

Unique, low edge, split shell design offers aseptic presentation protection that helps prevent accidental contamination of sets.

## **FEATURES AND BENEFITS**

- Advanced filter system includes puncture-proof filter protection. The system features interchangeable filter plates, offset baffle plates, easy-to-use filters and field replaceable main gaskets.
- TamperSTOP™, a patented integrated tamper-proof device and data card, combines three functions into one: process verification, data entry area and security.
- ACT2's space-saving design saves cost and shelf space, and requires only a one-inch clearance between the inner basket and the container.
- Color-coded silicone handle sleeves provide a secure, thermal-resistant grip.
- Full range of instrument protection and presentation accessories, including the Laparoscopic Rack, provide lightweight and compact solutions to your instrument handling needs.



**Teleflex®**  
MEDICAL



**J121740P**  
Medium Container

Pilling's ACT2 Sterilization Container Systems are crafted from fade- and scratch-resistant, clear hard-coated aluminum.



**B111130P**  
Small Basket

ACT2's space-saving design requires only a one-inch clearance between the inner basket and the container when fully assembled.



Identification Tag and Silicone Handle

Identification tags and silicone handles, available in a variety of colors to match sets, are offered free of charge with your Pilling ACT2 container purchase.



**o87403**  
UltraSoft™ Silicone Mat

Ultra Soft Silicone Mats are available for use with all basket sizes – large, medium, small and mini.

## CONTAINERS

### Large

12" Width x 23" Length\*

Perf. Top & Bottom	Solid Bottom	Depth
J122340P	K122340P	4"
J122350P	K122350P	5"
J122360P	K122360P	6"
J122370P	K122370P	7"

### Medium

12" Width x 17" Length

Perf. Top & Bottom	Solid Bottom	Depth
J121740P	K121740P	4"
J121750P	K121750P	5"
J121760P	K121760P	6"
J121770P	K121770P	7"

### Small

12" Width x 12" Length

Perf. Top & Bottom	Solid Bottom	Depth
J121240P	K121240P	4"
J121250P	K121250P	5"
J121260P	K121260P	6"
J121270P	K121270P	7"

### Mini

7 1/2" Width x 11" Length

Perf. Top & Bottom	Solid Bottom	Depth
Jo71140P	Ko71140P	4"
Jo71150P	Ko71150P	5"
Jo71160P	Ko71160P	6"

## CUSTOM SIZES

Custom sizes are available. Call with your requirements.



**o87362**

## STACKING BASKET ADAPTER

o87362 – Stacking Basket Adapters, set of four per package, are designed for stacking any basket on top of another, complete with locking nuts.

## BASKETS

### Large

11" Width x 22" Length

Basket No.	Depth
B112220P	2"
B112230P	3"
B112240P	4"
B112250P	5"
B112260P	6"

### Medium

11" Width x 16" Length

Basket No.	Depth
B111620P	2"
B111630P	3"
B111640P	4"
B111650P	5"
B111660P	6"

### Small

11" Width x 11" Length

Basket No.	Depth
B111120P	2"
B111130P	3"
B111140P	4"
B111150P	5"
B111160P	6"

### Mini

6 1/2" Width x 10" Length

Basket No.	Depth
Bo61020P	2"
Bo61030P	3"
Bo61040P	4"
Bo61050P	5"

\*All dimensions are nominal.

## PERFORATION TYPES

Two perforation types for baskets are available.

Perf. Type	Model No.
Standard	Use model # listed
Fine	Add suffix F, i.e., B112240FP

## IDENTIFICATION TAGS

Engraved aluminum tags are supplied **free of charge** with your Pilling ACT2 container purchase. Up to three lines of 15 characters each (including spaces) are provided. Colors available include: red, silver, blue, green, yellow, black, purple, orange and tan. *Note selection when ordering.* Also available, basket clip/tag A5021P.

## COLOR-CODED HANDLES

Extruded silicone handle sleeves provide a secure, thermal-resistant grip. Colors include: red, clear, blue, green, yellow, black, purple, orange and tan. *Note selection when ordering.*

## FILTERS

Disposable, 100 percent cellulose filters for both EtO and steam cycles are available, 1000 per case.

Size	Indicator	Model No.
9" x 6"	No	o83004
9 1/2" x 9 1/2"	No	o83001

## TAMPERSTOP™

A2004P – TamperSTOP indicator, all-in-one tamper-proof device/data card for *steam processing only*, 250 cycles per case.

## DATA CARDS

o82000 – Data Card, 250 per pack.

## TAMPER EVIDENT DEVICES

Color-coded Tamper Evident Devices are available with or without indicators.

	Qty	With Indicator	Without Indicator
Green	100	---	o81001
White	500	o81030	o81037
White	1000	o81031	o81034
Yellow	150	o81032	o81035
Yellow	500	o81033	o81036

## SILICONE MATS

Item No.	Description
o87402	6 1/2" x 10" UltraSoft mat
o87403	11" x 11" UltraSoft mat
o87419	11" x 16" UltraSoft mat
o87396	11" x 22" UltraSoft mat

## HOLDING PINS

Item No.	Description
o87349	3/4" Pins, set of 4
o87348	1 1/2" Pins, set of 6
o87350	1 1/2" Pin, 1 ea
o87341	2" Pin, 1 ea
o87351	3" Pin, 1 ea
o87352	3" Pins, set of 10

## DIVIDERS

Item No.	Description
o87353	(2) 5" x 3" Dividers w/ 10 ea 3" Pins
o87355	Corner Kit, 5" x 5" x 2 7/8" w/ 2 ea 2" Pins
o87356	Corner Kit, 6 1/2" x 6 1/2" x 3" w/2 ea 2" Pins
o87357	Corner Kit, 7 1/2" x 7 1/2" x 3" w/2 ea 2" Pins
o87359	(2) 6" x 3" Dividers w/ 4 ea 3" Pins
o87360	(2) 13" x 3" Dividers w/ 4 ea 3" Pins
o87361	(2) 15" x 3" Dividers w/ 4 ea 3" Pins
o87364	Corner Kit, 5" x 5" x 1 1/2" w/2 ea 1 1/2" Pins
o87365	Corner Kit, 3" x 3" x 1 1/2" w/2 ea 1 1/2" Pins
o87366	Corner Kit, 7" x 7" x 1 1/2" w/2 ea 1 1/2" Pins
o87367	(1) 1 1/2" x 6" Divider w/ 2 ea 1 1/2" Pins
o87368	(1) 1 1/2" x 10" Divider w/ 2 ea 1 1/2" Pins
o87369	(1) 1 1/2" x 15" Divider w/ 2 ea 1 1/2" Pins
o87371	(1) 2 1/4" x 6 1/2" Divider w/ 2 ea 2" Pins
o87372	(1) 2 1/4" x 11" Divider w/ 2 ea 2" Pins
o87373	(1) 2 1/4" x 13" Divider w/ 2 ea 2" Pins
o87374	(1) 2 1/4" x 17" Divider w/ 2 ea 2" Pins

## LAPAROSCOPIC RACK

The Pilling Laparoscopic Rack is a lightweight, space efficient solution to your laparoscopic handling needs. Get the maximum number of instruments into a compact and containerizable rack.

### FEATURES

Three level interlocking rack system is containerizable using Pilling's J122370P (perforated top and bottom) for sterilizing instruments without wrapping every time.

Top and middle deck hold any combination of 15 (DIN) or 17 (Standard) 5mm or 10mm ring handled or straight laparoscopic instruments.

Unique design allows for proper instrument sterilizing and drying with advanced open air construction for improved sterilization efficiency.

Constructed of durable aluminum and latex-free silicone materials.

Purchase each deck separately for easy organization of sets.



Pilling's Laparoscopic Rack is lightweight, space efficient and functional, providing a compact solution to your laparoscopic instrument handling needs.

## LAPAROSCOPIC RACK SIZES

### Standard

181120MCP	Three-tier Laparoscopic Rack, nominal outer dimensions: 11" W x 22" L x 6" H
181120MCTP	Top deck only, holds up to nine 5mm or 10mm diameter instruments
181120MCMP	Middle deck only, holds up to eight 5mm or 10mm diameter instruments
181120MCBP	Bottom deck only, holds miscellaneous instrumentation or accessories

Compatible: Complete set can fit into Pilling's ACT2™ Sterilization Container J122370P (perforated bottom) or K122370P (solid bottom).

### DIN

181020MNP	Complete three-piece set, nominal outer dimensions: 10" W x 21" L x 7" H
181020MNTP*	Top deck only, holds up to eight 5mm or 10mm diameter instruments
181020MNMP	Middle deck only, holds up to seven 5mm or 10mm diameter instruments
181020MNBP	Bottom deck only, holds miscellaneous instrumentation or accessories
181020MTMP	Top and middle decks
181020MMBP	Middle and bottom decks

Compatible: Complete set can fit into appropriately sized Aesculap® or comparable DIN-sized containers.

\*Top deck cannot be used alone with bottom deck. To use top and bottom deck you must use the middle deck.

\*Aesculap is a registered trademark of Aesculap Corporation.

## **PILLING**

Founded in 1814, the world's premier surgical instrument manufacturer is dedicated to designing and manufacturing hand-held surgical instruments of the highest caliber. Surgeons worldwide, in open and minimally invasive procedures, prefer Pilling instruments.

## **WARRANTY**

Lifetime Warranty against defects in workmanship or material on aluminum hard-coat, filter retention plate, baskets, accessories, handle assembly (excluding the actual handle), and posts. (These posts have been pull-tested to withstand up to 400 p.s.i. upward force without weld breakage occurring).

Three (3) Year Warranty or 200 cycles for inner and supplemental baffle plate gaskets.

Exclusions include:

- Misuse and abuse of the product
- Normal wear and tear of product
- Repairs completed at the hospital by other than factory-authorized personnel

## **ORDERING OPTIONS**

Order direct from your Pilling sales professional, or contact Pilling Customer Service:

**Phone: 800.523.6507**

**Fax: 800.332.2308**

International inquiries, contact:

**Phone: +919.361.3961**

**Email: [international@tfxsurgical.com](mailto:international@tfxsurgical.com)**

Or visit us online at:

**[www.pillingsurgical.com](http://www.pillingsurgical.com)**

Lit No. 990525

## **TELEFLEX**

Teleflex Medical is a division of Teleflex Incorporated (NYSE: TFX), a diversified global organization with annual revenues of approximately \$2 billion. Teleflex designs, manufactures and distributes quality-engineered products and services for the marine, industrial, automotive, aerospace, power generation and medical markets worldwide. The company employs more than 18,000 people worldwide who focus on providing innovative solutions for customers.